



Promise™ System for Proton Therapy

- Lightweight and Durable
- Comfortable Ergonomic Design
- Simple and Straightforward Mask Placement



Klarity Promise™ Baseplate
(Available for CT sim,
MR sim and treatment delivery)

The Klarity Promise™ Proton System features a lightweight baseplate designed for proton beam compatibility. The baseplate is equipped with comfortable, easily adjustable ergonomic hand grips.

The Promise System is compatible with Proton-K Masks, which feature our signature, non-stick Klarity Green® thermoplastic.

In addition to a standard carbon-fiber overlay, the Promise System is also available in MR compatible fiberglass as well as an insert-style for direct couch top treatments.





- R680-ACF
Promise™ System**
- 121 x 53 cm
 - 2.4 kg
 - Carbon Fiber
 - Equipped with ergonomic hand grips
 - Compatible with Karity Proton-K Masks




- R680-AGF
MR Promise™ System**
- 121 x 53 cm
 - Fiberglass
 - Equipped with ergonomic hand grips
 - Compatible with Karity Proton-K Masks

- R680-ACFP1
Promise™ System Couch Insert**
- 157 x 53 cm
 - Carbon Fiber
 - Attaches to couch top via pin-locks
 - Equipped with ergonomic hand grips
 - Compatible with Karity Proton-K Masks



**RG-PH-2436C
Green Proton-K Masks**
with 36% perforation
2.4mm

**RG-PH-3242C
Green Proton-K Head Mask**
with 42% perforation
3.2mm



**RG-PS-2436C
Green Proton-K Head & Shoulder Mask**
with 36% perforation
2.4mm

**RG-PS-3242C
Green Proton-K Head & Shoulder Mask**
with 42% perforation
3.2mm

**R550-PS
Klarity Cushion™ for
the Promise™ System**

